

## Pulse Lavage System

### System Components

#### PERFORMANCE SPECIFICATIONS:

##### PULSE RATE

- Variable to 1,700 pulses per minute by trigger and pressure adjustment

##### FLOW RATE

- Variable to 0.9 liters per minute

##### OPERATING PRESSURE

- 75 to 110 psi (5 - 7bar)

- Ergonomically designed and balanced, the L31 Pulse Lavage is an easy-to-use, pneumatic instrument that pumps fluids in distinct pulses at optimum pressure to remove blood and debris.
- The L31 Pulse Lavage features an easy-to-load pump cartridge to facilitate easy operation.
- A wide range of accessories including disposable tubing and nozzle sets designed for use during general purpose, hip and knee replacement surgeries are available.
- Unique groove in the Pulse Lavage handle keeps the delivery tube neatly tucked in and out of the way.
- Available with many popular hose connectors.

▲ Shown above with the 10.25"(25cm) femoral nozzle and suction attachment (Product No. L31-046).

Tava Surgical  
Pulse Lavage  
System

37



## Pulse Lavage System

### Product Number

### Hose Connection

L31-100	Compatible to a Tava Surgical Hose Connection
L31-200	Compatible to a Linvatec® Hall® Hose Connection
L31-300	Compatible to a 3M® Hose Connection
L31-400	Compatible to a Stryker® Hose Connection
L31-500	Compatible to a Synthes® Hose Connection
L31-600	Compatible to a Aesculap® Hose Connection

Tava Surgical  
Pulse Lavage  
System



▲  
Shown above with the knee arthroplasty set  
(Product No. L21-032).



# Pulse Lavage System


## System Accessories

<b>A90-000</b>	Air Hose compatible to a Tava Surgical/MicroAire® Connector (Length: 3m/10 feet)
<b>A90-000-15</b>	Air Hose compatible to a Tava Surgical/MicroAire® Connector (Length: 5m/15 feet)
<b>A90-107</b>	Air Hose (UK MK7 Style) compatible to a Tava Surgical/MicroAire® Connector (Length: 3m/10 feet)
<b>A90-107L</b>	Air Hose (UK MK7 Style) compatible to a Tava Surgical/MicroAire® Connector (Length: 5m/15 feet)
<b>A90-200</b>	Air Hose compatible to a Linvatec® Hall® Connector (Length: 3m/10 feet)
<b>A90-001</b>	Cleaning Plug/Cap compatible to a Tava Surgical/MicroAire® Connector
<b>A90-201</b>	Cleaning Plug/Cap compatible to a Linvatec® & Hall® Connector
<b>A90-501</b>	Cleaning Plug/Cap compatible to a A.O. Synthes® Connector
<b>A90-601</b>	Cleaning Plug/Cap compatible to a Aesculap® Connector
<b>M001</b>	Instrument Lube


Tava Surgical  
Pulse Lavage  
System

### Attachment Sets


<b>L21-030</b>	<b>L21-031</b>	<b>L21-032</b>
General Purpose Set	Hip Arthroplasty Set	Knee Arthroplasty Set



- Pump Cartridge
- 3m of tubing
- Bag spike connector
- 10cm shower head



- Pump Cartridge
- 3m of tubing
- Bag spike connector
- 10cm shower head acetabulum nozzle
- 25cm femoral shaft nozzle



- Pump Cartridge
- 3m of tubing
- Bag spike connector
- 10cm angled shower head



# Pulse Lavage System

## System Accessories



**L21-004**

10cm Showerhead Nozzle with Locking Collar



**L21-007**

22.5cm Femoral Shaft Nozzle with Locking Collar




**L21-005**

10cm Angled Showerhead Nozzle with Locking Collar




**L31-040**

Femoral Suction Attachment



**L31-042**

Suction Tubing




**L31-046**

Suction Tubing and Femoral Suction Attachment



**L21-010**

Splash Shield, 15cm Round



**L21-013**

Splash Shield, 12.5cm x 20cm Rectangular



**L21-017**

Cone (Cup) Splash Shield Ø14cm



**L21-020**

Cone (Cup) Splash Shield Ø14cm

